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1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 794-806.
2. Johnson AM, Rohlfes EM, Silverman LM. Proteins. In: Burtis CA, Ashwood ER. editors. Tietz textbook of clinical chemistry. 3rd ed. Philadelphia: W. B. Saunders Company; 1999. p. 507-12.
3. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996;34:517-20.

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BT 3500/1500/3000/2000/1000

| BT 3500/1500 | | |
|----------------------|-----------------------|-----------------|
| C4 | | |
| Primary Parameters | method | Sample Blank(A) |
| | Kind of process | Multi ponit |
| | filter | 340/- |
| | Reaction direction | increase |
| | REAGENT | |
| | Num of reagent | 2 |
| | Reagent1 | 180 |
| | Reagent2 | 30 |
| | SAMPLE | |
| | sample vol | 30 |
| | Pre dilution | 1 |
| | dilution | 5 |
| | TIME | |
| | Sample starter | inactive |
| | Delay time(sec) | 0 |
| | Reading time(sec) | 10 |
| Incubation time(sec) | 300/300 | |
| Check parameters | Reagent limit(M ABS) | 2000 |
| | Curve acceptance | 100 |
| | Test limit(conc) | 90 |
| | Initial ABS(M ABS) | 0 |
| | Final ABS(M ABS) | 3000 |
| | Max ABS delta (M ABS) | 2000 |
| | Chec prozone | inactive |
| | NORMAL RANGE | |
| | Male | 10 - 40 |
| | female | 10 - 40 |
| | Child | 10 - 40 |
| | Re-run hyperactive | inactive |
| | Re-run phatological | inactive |
| | Secondary | Unit serum |
| Num of needle washes | | 1/1 |
| Num of cuvette wash | | 2 |
| Additional wash | | inactive |
| Only needle | | |
| Only wash | | |
| Instrument factor | | 1 |
| Shift: | | 0 |
| Reagent blank | | everyrun |
| Decimals | | 0 |

| BT 3000/2000/1000 | | |
|----------------------|-----------------------|-----------------|
| C4 | | |
| Primary Parameters | method | Sample Blank(A) |
| | Kind of process | Multipoint |
| | filter | 340/- |
| | Reaction direction | increase |
| | REAGENT | |
| | Num of reagent | 2 |
| | Reagent1 | 300 |
| | Reagent2 | 50 |
| | SAMPLE | |
| | sample vol | 50 |
| | Pre dilution | 1 |
| | dilution | 5 |
| | TIME | |
| | Sample starter | inactive |
| | Delay time(sec) | 0 |
| | Reading time(sec) | 10 |
| Incubation time(sec) | 300/300 | |
| Check parameters | Reagent limit(M ABS) | 2000 |
| | Curve acceptance | 100 |
| | Test limit(conc) | 90 |
| | Initial ABS(M ABS) | 0 |
| | Final ABS(M ABS) | 3000 |
| | Max ABS delta (M ABS) | 2000 |
| | Chec prozone | inactive |
| | NORMAL RANGE | |
| | Male | 10 - 40 |
| | female | 10 - 40 |
| | Child | 10 - 40 |
| | Re-run hyperactive | inactive |
| | Re-run phatological | inactive |
| | Secondary | Unit serum |
| Num of needle washes | | 1/1 |
| Num of cuvette wash | | 2 |
| Additional wash | | inactive |
| Only needle | | |
| Only wash | | |
| Instrument factor | | 1 |
| Shift: | | 0 |
| Reagent blank | | everyrun |
| Decimals | | 0 |

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