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BT 3500/1500/3000/2000/1000

BT 3500/1500		
GGT		
Primary Parameters	method	kinetic
	Kind of process	linear
	filter	405/700
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	160
	Reagent2	40
	SAMPLE	
	sample vol	15
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
Reading time(sec)	120	
Incubation time(sec)	300/60	
Check parameters	Reagent limit(M ABS)	1000
	Curve acceptance	100
	Test limit(conc)	250
	Initial ABS(M ABS)	1100
	Final ABS(M ABS)	1600
	Max ABS delta (M ABS)	500
	Chec prozone	inactive
	NORMAL RANGE	
	Male	0 - 49
	female	0 - 32
	Child	-
	Re-run hyperactive	inactive
	Re-run phatological	inactive
Secondary	Unit serum	U/L
	Num of needle washes	1/1
	Num of cuvette wash	2
	Additional wash	inactive
	Only neadle	
	Only wash	
	Instrument factor	1
	Shift:	0
	Reagent blank	everyday
	Decimals	0

BT 3000/2000/1000		
GGT		
Primary Parameters	method	kinetic
	Kind of process	linear
	filter	405/700
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	240
	Reagent2	60
	SAMPLE	
	sample vol	24
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
Reading time(sec)	120	
Incubation time(sec)	300/60	
Check parameters	Reagent limit(M ABS)	1000
	Curve acceptance	100
	Test limit(conc)	250
	Initial ABS(M ABS)	1100
	Final ABS(M ABS)	1600
	Max ABS delta (M ABS)	500
	Chec prozone	inactive
	NORMAL RANGE	
	Male	0 - 49
	female	0 - 32
	Child	-
	Re-run hyperactive	inactive
	Re-run phatological	inactive
Secondary	Unit serum	U/L
	Num of needle washes	1/1
	Num of cuvette wash	2
	Additional wash	inactive
	Only neadle	
	Only wash	
	Instrument factor	1
	Shift:	0
	Reagent blank	everyday
	Decimals	0

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