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BT 3500/1500/3000/2000/1000

BT 3500/1500		
Iron		
Primary Parameters	method	Sample Blank(A-B)
	Kind of process	linear
	filter	578/-
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	180
	Reagent2	45
	SAMPLE	
	sample vol	40
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
Reading time(sec)	10	
Incubation time(sec)	300/300	
Check parameters	Reagent limit(M ABS)	2000
	Curve acceptance	100
	Test limit(conc)	250
	Initial ABS(M ABS)	0
	Final ABS(M ABS)	3000
	Max ABS delta (M ABS)	2000
	Chec prozone	inactive
	NORMAL RANGE	
	Male	35 - 168
	female	37 - 165
	Child	22 - 201
	Re-run hyperactive	inactive
Re-run phatological	inactive	
Secondary	Unit serum	Mic g/dl
	Num of needle washes	2/2
	Num of cuvette wash	3
	Additional wash (Acid)	active
	Only neadle	inactive
	Only one wash	active
	Instrument factor	1
	Shift:	0
	Reagent blank	everyrun
	Decimals	0

BT 3000/2000/1000		
Iron		
Primary Parameters	Method	Sample Blank(A-B)
	Kind of process	linear
	Filter	578/-
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	300
	Reagent2	75
	SAMPLE	
	sample vol	60
	Pre dilution	1
	Dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
Reading time(sec)	10	
Incubation time(sec)	300/300	
Check parameters	Reagent limit(M ABS)	1000
	Curve acceptance	100
	Test limit(conc)	250
	Initial ABS(M ABS)	0
	Final ABS(M ABS)	3000
	Max ABS elta (M ABS)	2000
	Chec prozone	inactive
	NORMAL RANGE	
	Male	35 - 168
	female	37 - 165
	Child	22 - 201
	Re-run hype ractive	inactive
Re-run phatological	inactive	
Secondary	Unit serum	Mic g/dl
	Num of needle washes	2/2
	Num of cuvette wash	3
	Additional wash (Acid)	active
	Only neadle	inactive
	Only one wash	active
	Instrument factor	1
	Shift:	0
	Reagent blank	everyrun
	Decimals	0

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