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BT 3500/1500/3000/2000/1000

BT 3500/1500		
RF		
Primary Parameters	method	Sample Blank(A)
	Kind of process	linear
	filter	340/-
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	180
	Reagent2	18
	SAMPLE	
	sample vol	24
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
	Reading time(sec)	10
	Incubation time(sec)	300/300
Check parameters	Reagent limit(M ABS)	2000
	Curve acceptance	100
	Test limit(conc)	40
	Initial ABS(M ABS)	3000
	Final ABS(M ABS)	2000
	Max ABS delta (M ABS)	300
	Chec prozone	inactive
	NORMAL RANGE	
	Male	0 - 15
	female	0 - 15
	Child	0 - 15
	Re-run hyperactive	inactive
	Re-run phatological	inactive
Secondary	Unit serum	IU/ml
	Num of needle washes	1/1
	Num of cuvette wash	2
	Additional wash	inactive
	Only neadle	
	Only wash	
	Instrument factor	1
	Shift:	0
	Reagent blank	everyday
	Decimals	1

BT 3000/2000/1000		
RF		
Primary Parameters	method	Sample Blank(A)
	Kind of process	linear
	filter	340/-
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	300
	Reagent2	30
	SAMPLE	
	sample vol	40
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
	Reading time(sec)	10
	Incubation time(sec)	300/300
Check parameters	Reagent limit(M ABS)	2000
	Curve acceptance	100
	Test limit(conc)	40
	Initial ABS(M ABS)	3000
	Final ABS(M ABS)	2000
	Max ABS delta (M ABS)	300
	Chec prozone	inactive
	NORMAL RANGE	
	Male	0 - 15
	female	0 - 15
	Child	0 - 15
	Re-run hyperactive	inactive
	Re-run phatological	inactive
Secondary	Unit serum	IU/ml
	Num of needle washes	1/1
	Num of cuvette wash	2
	Additional wash	inactive
	Only neadle	
	Only wash	
	Instrument factor	1
	Shift:	0
	Reagent blank	everyday
	Decimals	1

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