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Range	
Application Code	060
Report Name	CK-MB
Data Mode	On Board
Control Interval	1000
Instrument Factor (Y=aX+b)	a= 1.0 b= 0.0
Technical Limit	0
Repeat Limit	0
Expected Value	

Others	
<Standard>	(1) (2) (3) (4) (5) (6)
Calib. Code	501
Concentration	0
Position	*
Sample Volume	8
Diluent S.Vol	0
Diluent Volume	0

#) Data entry by the user
*) Enter calibration or standard value

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Imidazole	PH 6.7	100 mmol/l
Creatine phosphate		30 mmol/l
Glucose		20 mmol/l
N-Acetyl cysteine	(NAC)	20 mmol/l
Magnesium acetate		10 mmol/l
EDTA-Na ₂		2 mmol/l
ADP		2 mmol/l
NADP		2 mmol/l
AMP		5 mmol/l
Diadenosine pentaphosphate		10 µmol/l
Glucose-6-phosphate dehydrogenase	(G6P-DH)	≥ 1.5 KU/L
Hexokinase	(HK)	≥ 2.5 KU/L
Monoclonal antibodies (sheep) against human CK-M;inhibiting capacity		≥2000 U/L

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- Stein W. Creatine kinase (total activity), creatine kinase isoenzymes and variants. In: Thomas L, ed. Clinical laboratory diagnostics. Frankfurt: TH-Books Verlagsgesellschaft;1998.p.71-80.
- Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.
- Würzburg U, Hennrich N, Orth HD, Lang H. Quantitative determination of creatine kinase isoenzyme catalytic concentrations in serum using immunological methods. J Clin Chem Clin Biochem 1977;15:131-7.
- Recommendations of the German Society for Clinical Chemistry. Standardization of methods for the estimation of enzyme activities in biological fluids: Standard method for the determination of creatine kinase activity. J Clin Chem Clin Biochem 1977;15:255-60.

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