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100 mmol/l
PH 7.5
TRIS buffer
Polyethyleneglycol , detergents and stabilizers

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100 mmol/l
PH 8.0
TRIS buffer
Anti-human CRP antibodies (goat) with stabilizers

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1. Thompson D, Milford-Ward A, Whicher JT. The value of acute phase protein measurements in clinical practice. Ann Clin Biochem 1992;29:123-31.
2. Gabay C, Kushner I. Acute-phase proteins and other systemic responses to inflammation. N Engl J Med 1999;340:448-54.
3. Hansson LO, Lindquist L. C-reactive protein: its role in the diagnosis and follow-up of infectious diseases. Curr Opin Infect Diseases 1997; 10:196-201.
4. Sipe JD. Acute-phase proteins in osteoarthritis. Semin Arthritis Rheum 1995;25:75-86.
5. Ridker PM, Cushman M, Stampfer MJ, Tracy RP, Hennekens CH. Inflammation, aspirin and the risk of cardiovascular disease in apparently healthy men. N Engl J Med 1997;336:973-9.
6. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Biennu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). Eur J Clin Chem Clin Biochem 1996;34:517-20.

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BT 3500/1500/3000/2000/1000

BT 3500/1500		
CRP		
Primary Parameters	method	Sample Blank(A)
	Kind of process	linear
	filter	340/-
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	180
	Reagent2	18
	SAMPLE	
	sample vol	12
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
	Reading time(sec)	10
	Incubation time(sec)	300/300
Check parameters	Reagent limit(M ABS)	2000
	Curve acceptance	100
	Test limit(conc)	30
	Initial ABS(M ABS)	0
	Final ABS(M ABS)	3000
	Max ABS delta (M ABS)	2000
	Chec prozone	inactive
	NORMAL RANGE	
	Male	0 - 10
	female	0 - 10
	Child	-
	Re-run hyperactive	inactive
	Re-run phatological	inactive
Secondary	Unit serum	mg/L
	Num of needle washes	1/1
	Num of cuvette wash	2
	Additional wash	inactive
	Only neadle	
	Only wash	
	Instrument factor	1
	Shift:	0
	Reagent blank	everyrun
	Decimals	1

BT 3000/2000/1000		
CRP		
Primary Parameters	method	Sample Blank(A)
	Kind of process	linear
	filter	340/-
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	300
	Reagent2	30
	SAMPLE	
	sample vol	20
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
	Reading time(sec)	10
	Incubation time(sec)	300/300
Check parameters	Reagent limit(M ABS)	2000
	Curve acceptance	100
	Test limit(conc)	30
	Initial ABS(M ABS)	0
	Final ABS(M ABS)	3000
	Max ABS delta (M ABS)	2000
	Chec prozone	inactive
	NORMAL RANGE	
	Male	0 - 10
	female	0 - 10
	Child	-
	Re-run hyperactive	inactive
	Re-run phatological	inactive
Secondary	Unit serum	mg/L
	Num of needle washes	1/1
	Num of cuvette wash	2
	Additional wash	inactive
	Only neadle	
	Only wash	
	Instrument factor	1
	Shift:	0
	Reagent blank	everyrun
	Decimals	1

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