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1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH- Books Verlagsgesellschaft; 1998. p. 142-48.
2. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Titez Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B. Saunders Company; 1999. p. 790-6.
3. Hoelzel W, Weykamp C et al. IFCC Reference System for measurement of Hemoglobin A1c in Human Blood and the National Standardization Schemes in the United States, Japan, and Sweden: A Method comparison Study. Clin Chem 2004; 50:1:166-74.
4. Midema K. Standarization of HbA1C and Optimal Range of Monitoring. Scand J clin Lab Invest 2005/65 (suppl 240):61-72

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BT 3500/1500/3000/2000/1000

BT 3000/2000/1000		
HBA1C		
Primary Parameters	method	End point
	Kind of process	Multi ponit
	filter	630/-
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	190
	Reagent2	63(Dilutio solution32)
	SAMPLE	
	sample vol	5
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
	Reading time(sec)	10
	Incubation time(sec)	300/300
Check parameters	Reagent limit(M ABS)	2000
	Curve acceptance	100
	Test limit(conc)	23
	Initial ABS(M ABS)	0
	Final ABS(M ABS)	3000
	Max ABS delta (M ABS)	2000
	Chec prozone	inactive
	NORMAL RANGE	
	Male	4 - 6
	female	4 - 6
	Child	4 - 6
	Re-run hyperactive	inactive
	Re-run phatological	inactive
Secondary	Unit serum	%
	Num of needle washes	2/2
	Num of cuvette wash	2
	Additional wash	active
	Only neadle	active
	Only wash	active
	Instrument factor	1
	Shift:	0
	Reagent blank	everyday
	Decimals	2

BT 3500/1500		
HBA1C		
Primary Parameters	method	End point
	Kind of process	Multi ponit
	filter	630/-
	Reaction direction	increase
	REAGENT	
	Num of reagent	2
	Reagent1	120
	Reagent2	40(Dilution solution 20)
	SAMPLE	
	sample vol	4
	Pre dilution	1
	dilution	5
	TIME	
	Sample starter	inactive
	Delay time(sec)	0
	Reading time(sec)	10
	Incubation time(sec)	300/300
Check parameters	Reagent limit(M ABS)	2000
	Curve acceptance	100
	Test limit(conc)	23
	Initial ABS(M ABS)	0
	Final ABS(M ABS)	3000
	Max ABS delta (M ABS)	2000
	Chec prozone	inactive
	NORMAL RANGE	
	Male	4 - 6
	female	4 - 6
	Child	4 - 6
	Re-run hyperactive	inactive
	Re-run phatological	inactive
Secondary	Unit serum	%
	Num of needle washes	2/2
	Num of cuvette wash	2
	Additional wash	active
	Only neadle	active
	Only wash	active
	Instrument factor	1
	Shift:	0
	Reagent blank	everyday
	Decimals	2

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