

GPT/ALT		
primary parameters	Method	Kinetic
	Kind of Process	linear
	Filter	340 / 405
	Reaction Direction	Decrease
	reagent	
	Num of Reagent	2
	Reagent 1	240
	Reagent 2	60
	sample	
	Sample Vol	30
	Pre Dilution	1
	Dilution	5
	time	
	Sample Starter	Active
	Delay Time (sec)	0/30
	Reading Time (sec)	180
	Incubation Time (sec)	60/60
check parameters	Reagent Limit (M ABS)	0
	Curve Acceptance (%)	100
	Test Limit (Conc)	250
	Initial ABS (M ABS)	3000
	Final ABS (M ABS)	500
	Max ABS Delta (M ABS)	200
	Check Prozone	Inactive
	normal range	
	Male	2 - 38
	Female	2 - 38
	Child	2 - 38
	Re - Run Hyperactive	Active
	Re - Run Phatological	INACTIVE
secondary parameters	Unit Serum	U / L
	Num OF Needle Washes	1/1
	Num OF Cuvette Washes	1
	Aditonal Wash	Inactive
	Only neadle	
	Only Wash	
	Instrument Factor	0
	Shift :	0
	Reagent Blank	EVERY DAY
	Decimals	0