

Triglycerides		
primary parameters	Method	End Point
	Kind of Process	linear
	Filter	510 / 700
	Reaction Direction	INCREASE
	reagent	
	Num of Reagent	1
	Reagent 1	300
	Reagent 2	
	sample	
	Sample Vol	3
	Pre Dilution	1
	Dilution	5
	time	
	Sample Starter	Inactive
	Delay Time (sec)	0
	Reading Time (sec)	10
	Incubation Time (sec)	300
check parameters	Reagent Limit (M ABS)	2000
	Curve Acceptance (%)	100
	Test Limit (Conc)	1000
	Initial ABS (M ABS)	—
	Final ABS (M ABS)	—
	Max ABS Delta (M ABS)	600
	Check Prozone	Inactive
	normal range	
	Male	40 - 160
	Female	40 - 160
	Child	40 - 160
	Re - Run Hyperactive	Active
	Re - Run Phatological	INACTIVE
secondary parameters	Unit Serum	MG / DL
	Num OF Needle Washes	1
	Num OF Cuvette Washes	1
	Aditonal Wash	Inactive
	Only neadle	
	Only Wash	
	Instrument Factor	0
	Shift :	0
	Reagent Blank	EVERY DAY
	Decimals	0